

Ray's Fluid Thioglycollate Method*

A Standard Practical Guide for
Oyster Sentinel Participants

*as demonstrated by Sammy M. Ray

Preparing Dermo tubes:

1. Add 20 gm NaCl^a to 1L of deionized (DI) water.
2. Add 29.0 gm of thioglycollate^b to the water, heat on low temperature and mix to dissolve.
3. Dispense 10 ml of medium into screw cap culture tubes. (Have caps on loosely for autoclaving.)
4. Autoclave for 15 min.
5. Allow the tubes to cool then tighten the caps.
6. Store the tubes in the dark at room temperature until needed.

Suggested Vendors:

- a. Sigma 55886-500G, 500 grams of sodium chloride
- b. Fluka 70157 – 500G, 500 grams of Thioglycollate Broth



Preparing antibiotics:

Stock Nystatin^c (=Mycostatin)


1. Add 9ml of DI water to the 5 million unit vial and shake well.
2. Add 2.5 ml of the re-hydrated nystatin into each of 4 vials. Date and label as Nystatin Stock 1, 2, 3 and 4.
3. Freeze until needed. (Good for at least a year.)

Chloromycetin/Nystatin Working Solution:

1. Add 4.5 ml of DI water to a 1 gm vial of Chloromycetin^d (=Chloramphenicol) and shake well.
2. Add the re-hydrated Chloromycetin to the Nystatin Stock vial.
3. Add 17.5 of DI water to the Chloromycetin/Nystatin mix.
4. Date and label the vial as Chlor/Nystatin Working Solution.
5. Refrigerate.

Suggested Vendors:

- c. Sigma N6261 - 5MU, 5 Million unit vial of Nystatin
- d. Sigma C3738 – 5G, 5 grams of Chloramphenicol Succinate Sodium Salt



Inoculation of tubes:

1. Shake the Chloromycetin/Nystatin Working Solution well.
2. Add 0.05 ml of the Chloromycetin/Nystatin Working Solution to each Dermo tube and mix by inverting the tube.
3. Place oyster ~ 5mm² piece of anterior mantle tissue into Dermo tube. (Be sure the tissue is in the fluid!)
4. Store in the dark at room temperature for about a week.



Preparation of Lugol's working solution and staining of samples:

1. Prepare 50 ml of working solution by adding 40 ml of distilled or deionized water to 10ml of 1N Iodine^e stock solution.
2. With an inoculating needle, carefully remove the oyster tissue from the tube and place it on a glass slide.
3. Add 2-3 drops of the Lugol's working solution to the tissue.
4. Spread and macerate the tissue with a blunt probe to get a thin, well-stained preparation.
5. Cover the tissue with a cover slip. Push on the cover slip to get a flattened preparation.
6. Remove excess Lugol's with absorbent paper.



Suggested Vendors:

e. Fisher S178 – 500, 500 ml of 1N Iodine Solution

Reading of samples:

1. Determine level of parasitism using the Mackin (1962) 0-5 scale, as modified by Craig et al. (1989). Scan the slide under low magnification (e.g., 40x) then switch to 100x magnification if necessary.
2. Use the following photomicrographs as standards as needed. Note that the field of view shown in the following slides is intended to represent the entire sample.

References:

Craig, A., E. N. Powell, R. R. Fay & J. M. Brooks. 1989. Distribution of *Perkinsus marinus* in gulf coast oyster populations. *Estuaries*. 12:82-91.

Mackin, J. G. 1961. Oyster disease caused by *Dermocystidium marinum* and other microorganisms in Louisiana. *Publ. Inst. Mar. Sci.* 7:132-299.



Dermo Code:


0.00

Description:

No hyphospores present

Comment:

Note the numerous Brown
Cells , which are NOT Dermo



Dermo Code:


0.33

Description:

1-10 hyphospores

Comment:

Note the single Dermo cell
among the many Brown Cells



Dermo Code:


0.67

Description:

11-74 hypospores

Comment:

Note the dark-stained
spherical Dermo cells
among the Brown Cells



Dermo Code:


1.00

Description:

75-125 hypospores

Comment:

Note the many Dermo cells
among the many Brown Cells



Dermo Code:


1.33

Description:

>125 hyphospores but much less than 25% of tissue is hyphospores. Do not count beyond 125 hyphospores

Comment:

Note the slight variation in hyphospore size and staining intensity



Dermo Code:


1.67

Description:

<25% of tissue is hyphospores

Comment:

Starting with code 1.67,
dermo codes are based on
percent coverage not numbers



Dermo Code:

2.00


Description:

25% of tissue is hyphae

Comment:

It is difficult to judge an exact percent coverage such as this. It helps to judge the sample in relationship to the codes below and above.

Note also that dermo is patchy and that the percent coverage must be judged for the entire sample.



Dermo Code:


2.33

Description:

>25% but much less than 50%
of tissue is hyphospores

Comment:

Note the patchy distribution of
the hyphospores



Derma Code:

2.67


Description:

>25% but <50% of tissue is
hyphospores

Comment:

Some areas are >50% covered,
whereas other areas are <25%
covered.

A judgment is made on the
average percent coverage of
the whole sample.



Dermo Code:

3.00


Description:

50% of tissue is hyphospores

Comment:

It is difficult to judge an exact percent coverage such as this. It helps to judge the sample in relationship to the codes below and above.

Note also the variation in hyphospore size which affects percent coverage.



Derma Code:

3.33


Description:

>50% but *much less* than 75%
of tissue is hyphospores

Comment:

Note the patchiness of
hyphospore distribution

Determination of disease code
must be made for the entire
sample which in this case
includes dense and sparse
patches



Derma Code:

3.67


Description:

>50% but <75% of tissue is
hyphospores

Comment:

Note the patchiness of
hyphospore distribution

Determination of disease code
must be made for the entire
sample, which in this example
includes dense and sparse
patches



Dermo Code:

4.00

Description:


75% of tissue is hyphospores

Comment:

Note the patchiness of hyphospore distribution.

It is difficult to judge an exact percent coverage such as this.

It helps to judge the sample in relationship to the codes below and above.



Dermo Code:


4.33

Description:

>75% but *much less* than 100%
of tissue is hyphospores

Comment:

This sample has a fairly even
distribution of hyphospores,
which however are of varying
size and staining intensity



Dermo Code:


4.67

Description:

>75% tissue is hyphospores
but some oyster tissue is still
visible

Comment:

Note the heavy coverage by
dermo but with oyster tissue
still visible



Dermo Code:


5.00

Description:

Nearly 100% of tissue is
hypnospores

Comment:

Note that the oyster tissue is
nearly totally obscured by
hypnospores



Data entry:

Go to

www.oystersentinel.cs.uno.edu

Enter data

Check data entry

Submit data

The web site will display:

Percent Infection as

$(\text{no. infected}/\text{no. assayed}) \times 100$

Infection Intensity as

$\text{sum of disease codes}/\text{no. infected}$

Weighted Prevalence as

$\text{sum of disease codes}/\text{no. assayed}$



Check the web site to see if
data were properly posted